



Experience the Trace Quality

*Happier Living Everyday*

SINGLE CHANNEL ECG RECORDER CARDIART 6108T

**Single Channel ECG Recording**

Enables the recording of all leads sequentially

**Ergonomic Design**

The compact, light weight and portable design ensures ease of use

**Automatic and Manual Recording Modes**

Manual mode enables the acquisition and printing of ECG recording in real time. Simultaneous acquisition is possible with automatic mode

**Last Memorised ECG Recording**

Ensures the storage of last recording in the memory

**Thermal Printing**

Allows the printing of ECG on 50 mm paper using a high-resolution thermal print

**Power Supply**

Comes with dual power supply – mains and rechargeable internal battery

**Built-in Rechargeable Battery**

Enables mains independent operations and provides battery use after charging

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



A BPL Medical Technologies
Magazine



Follow Us on

