



Experience
the Accuracy

MULTI PARAMETER MONITOR
EXCELLO



Happier Living Everyday



The device provides arrhythmia detection, ST & HRV analysis with 3-lead/5-lead ECG monitoring and comprehensive trend data review



High Resolution Display

High resolution 10.4" TFT colour display.



Standard Configuration

Standard configuration of ECG, RESP, Nellcor SpO₂ / Masimo SpO₂ / Digital SpO₂, NIBP, TEMP



Calculation

The added functions of the device are Drug dose calculation, HRV analysis & Hemodynamic calculation



Connectivity to CNS

Connectivity to optional BPL Medius+ Central Nursing Station Software (Wired or optional wireless)



168 Hour Trend

Provides for maximum of 168 hours of trend review



ECG Monitoring

5 Lead ECG / 3 Lead ECG (optional)



Waveform Display

Maximum 10 waveforms display



Arrhythmia and ST Analysis

The device provides for arrhythmia detection and ST analysis



Optional Parameters & Configuration

2-IBP, ETCO₂ (Mainstream / Side stream), Multi-Gas/ AGM, Noninvasive Cardiac Output, HL7 Output and Thermal Recorder