

Continuous Patient Care with **Efficient Monitoring**



MULTIPARAMETER PATIENT MONITOR SYSTEM
BPL VIVID VUE M12

Follow us on @bplmedtec h




www.bplmedicaltechnologies.com



MULTIPARAMETER PATIENT MONITOR SYSTEM

BPL VIVID VUE M12

With advanced monitoring capabilities and intuitive design, our patient monitor supports healthcare professionals in delivering exceptional patient care.



*Optional Connectivity to BPL Medius+ Central Nursing Station software



OxyCRG Screen

OxyCRG graphs HR, SPO₂ and Respiratory waveform continuously for easy comparison of breathing cessations & HR decelerations



Review

Trend and event can be reviewed with associated waveform and table values. Provision to print the same is available



QT View

QT View displays the current QT parameter values and waveforms as well as baseline QT parameter values and waveforms. Options to export QT View report is also available.



ST Analysis window

Along with standard automated ST monitoring, Vivid Vue M12 also has manual ST window where the user can set E and J isoelectric positions for ST analysis and also provide ST bar graph along with ST values


CO₂ Analyzer



CO₂ Analyzer Sidestream

Model: Masimo ISA CO₂

Measures: EtCO₂, FiCO₂, RR



BPL Capnoview+ Sidestream

Model: BPL Capnoview+


Measures: EtCO₂, FiCO₂, RR



Mainstream CO₂ Analyzer

Model: Masimo IRMA CO₂

Measures: EtCO₂, FiCO₂, RR




BPL Capnoview+ Mainstream

Model: BPL Capnoview+


Measures: EtCO₂, FiCO₂, RR

BPL iCARE OT: A digital solution that simplifies operating theatre management by automating vital signs, documenting events, and enabling remote monitoring.

BPL iCare OT
OT Charting Software Application



A Global Medical Technology Company



- BPL MedTech HQ (Bengaluru)
- BPL MedTech Factory (Bengaluru)
- BPL MedTech Factory (Palakkad)

Service Powered by

BPL SmartUKeep®
Service

Service Helpline
1800-425-2355
9AM - 8PM (Monday - Saturday)

ISO 13485:2016 CERTIFIED COMPANY
ISO 9001:2015 CERTIFIED COMPANY

BPL MEDICAL TECHNOLOGIES PRIVATE LIMITED

Regd. Office: 11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076, India.
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in
CIN: U33110KA2012PTC067282

EPR Compliance
BPL Medical Technologies Pvt. Ltd.
Thickness of plastic packaging is 50 microns and above
EPR Registration Number: IM-01-000-12-AAFCB3158E-23



Good at Heart®
A BPL Medical Technologies Initiative

Wake up to Wellness®
A BPL Medical Technologies Initiative

