



# Experience the Precision



## MULTI PARAMETER MONITOR

# EliteView EV10

*Happier Living Everyday*



### High Resolution Touch Screen Display

10.4" color LCD Touch Screen with high pixel resolution



### Standard Configuration

Increased clinical confidence in measuring the basic vital parameters - ECG, SPO2 (Masimo / Nellcor / Digital), NiBP, RESP, TEMP



### Optional Parameters & Configuration

2-IBP, ETCO2 (Main Stream / Side Stream), 12 Lead ECG and Thermal Recorder



### Calculation Interface

Drug Dose Calculation, Hemodynamic, Ventilation, Oxygenation and Renal Calculation provides users with more information regarding the status of the patient



### Enhanced Connectivity Solution

Central nursing station, USB patient data output, External printer, HL7, FTP server upload, VGA for external display



### Intelligent Alarms

Audio & Visual Alarms with 3 different levels of Priority: High, Medium & Low can help clinicians easily understand the criticality of patient condition



### High performance Li-ion Battery

With low power consumption, the battery back-up time is 2 hrs. (4hrs back-up can be provided as an option)



### Comprehensive Data Storage

160 Hours of Trend Data, 2000 groups of NiBP, 200 Alarm Review, 200 Arrhythmia recall, 48 hrs Full disclosure review; provides users for detailed review of patient status