

# IBIS



## iCore

Infant Warmer

# MEDICAL

Every newborn deserves the safest and most nurturing start to life. The iCore 502 infant warmer from Ibis Medical provides exceptional supplementary heat to the baby, with precise monitoring and utmost safety.

Our warmers provide infants with precise warmth they need for a safe and comfortable start to life. Engineered for both performance and ease of use, iCore 502 provides a unique blend of performance and safety.

### Safe

- iCore 502 comes with a unique first in class Sensor Displacement Detection System. If by any chance, the temperature sensor is displaced from patient body, iCore 502 will automatically regulate the temperature to 50% heater output & sound an alarm. Thus ensuring that the baby is safe at all times.
- iCore 502 infant warmers provide accurate audio visual alarms for easy problem identification.
- Baby safety is ensured with our prompt heater cut-off if the baby skin temperature exceeds 38° C.

### Feature Rich

iCore 502 comes with a host of useful features. These include

- A pre-warm mode in addition to Servo and Manual modes.
- APGAR timer.
- Timer with alarm which can be set between 30 minutes to 5 hours in half an hour intervals.
- Centigrade to Fahrenheit toggle button.
- All important audio visual alarms which include
  - High and low temperature alarms.
  - Power failure, probe failure, heater failure alarms.
  - High temperature cutoff alarm if the baby temperature exceeds 38° C.
  - Check Patient alarm, which sounds every 15 minutes (in manual mode)

### Engineered for Efficiency

iCore 502 infant warmers are well designed and engineered to do what matters most to you. Which is, to provide the optimum & uniform heating across the baby mattress, allowing greater thermal comfort to the baby.

- Our 600W quartz infrared heater coupled with parabolic design SS reflector provides fast and effective heating.
- Our super responsive thermistor sensors are specially designed only for our infant warmers. They are very fast and very accurate, which makes our warmers very sensitive to quick temperature changes on baby body.

## 502

Integrated Bed & Drawer